

Glow User Manual

Getting Started:

1. Allow Camera Permission when prompted.
2. Hold your phone at eye level and align your face within the white oval.
3. The oval will turn green when your face is detected.

Taking a Check-In:

1. Tap "Start 15s Check-In".
2. The AI voice assistant will guide you.
3. The screen will flash bright white for 5 seconds. Do not close your eyes.
4. After 15 seconds, your mental fatigue result (Rested, Fatigued, or Stressed) will appear.

Recommendations:

1. After saving your result, you will be taken to the Recommendations screen.
2. Perform suggested activities like Breathing, Mental Reset, or Stretching.
3. Click "Finish & Save Session" to record your activities.

History & Export:

1. Tap "  View History" on the home screen.
2. Tap any past check-in to add a Private Note.
3. Tap "  Export History as PDF" to generate a report for your doctor/therapist.